



**MICROGRAPH**  
TECHNOLOGIES

**Millipore Testing and Particle Analysis Service  
as per standard ISO 16232/VDA 19/and  
many more**



### **MICROSCOPE WITH MOTORIZED STAGE**

Micrograph Technologies can provide either the equipment or the testing service you require for the determination of the technical cleanliness of your components in accordance with inspection standards ISO16232




Contamination removal from components used in automotive, hydraulics and production systems is becoming increasingly important as suppliers, producers and end users look for better reliability.

The microscope analysis identifies the type & quantity of particles in different size classes including longest particle and fibres according to the applicable test standards applied.

The microscope system can separate particles into metallic and non-metallic to give a distribution of metallic and non-metallic particles.

### FLUID CONTAMINATION ANALYSIS KIT



vacuum Pressure Pump, Oil Free, Neoprene

Diaphragm coated with Teflon

HYDROSOL, SS Filter Holder, 47mm

Vacuum Hose, Silicone, 4.5" Long

Vacuum Filter Flask, Glass, 1 Ltr

- ▶ SS Forceps, Smooth Tip
- ▶ Solvent Filter Dispenser Bottle, 25mm
- ▶ Pressure Vessel, SS 316, 5 Ltr
- ▶ Filter Solvent Dispensing Gun, 25mm
- ▶ Vacuum Trap Assembly
- ▶ Petri dish.



## MICROSCOPE



### Stereo Zoom Microscope

**Magnification:** 7x to 45x

**Zooming Ratio:** 6.4:1

**Focussing Stroke:** 50mm

**Objective:** 0.7x to 4

**Eyepiece:** 10x (Pair)

**Illumination:** White Ring Light LED

### Trinocular Upright Metallurgical Microscope

**Magnification:** 50x to 1000x

**Trinocular Observation**

**Head Eyepiece (Anti Fungus)**

**Wide Field 10x (paired) F.O.V. 22mm**

**Illuminator:**

The EPI Illuminator has 6V-20W

halogen lamp with variable brightness control



## MOTORISED STAGE

**Motorized Stage:** Stage Size - 150mm X 150mm


**XY Movement:** 60mm X 60 mm (Motorized)

**Resolution:** 05 micron or better

**Stage Control:** Through software.



## PARTICLE ANALYSIS SOFTWARE



- ▶ Pre-configured for ISO 16232.
- ▶ Calculation of Component-Cleanliness-Code (CCC) according to ISO 16232
- ▶ Software controlled Stage movements.
- ▶ Auto image stitching and analysis.
- ▶ Automatic scanning and capturing the entire field of filter paper.
- ▶ Software generate a single image of entire filter paper (size 47mm).
- ▶ So asto analyse and deliver the results for their true sizes.
- ▶ Separation of particles as Metallic, Non-metallic and fibre.
- ▶ Measurement of all the particles for their total number, size, perimeter, area.
- ▶ Distribution Metallic and Non-metallic particles as per their sizes in user defined size classes.
- ▶ Automatic detection of particle types Reflective, Non-Reflective, Reflective- Fibrous and Non-Reflective-Fibrous.
- ▶ Analysis of particle parameters for each particle (max particle size, min particle size, length, width, area fibre length, and many more).




# PARTICLE ANALYSIS REPORT

<b>Customer Name:</b>	<b>Date:</b> 11/06/2020
<b>Part Name:</b> Piston Brake	<b>Drp Rev No:</b>
<b>Part No:</b> 10123003	<b>Drp Rev Date:</b>
<b>Ref Insp. Std:</b>	<b>Prod. Batch ID:</b>
<b>Test Method:</b> ISO 16232 / ISO 4407	<b>Microscope ID:</b>
<b>Paper Size:</b> 47 mm dia	<b>Magnification:</b> 50x
<b>Calibration Factor:</b> 2.4 micron / pixel	<b>No of Parts for Test:</b>
<b>No of Fields:</b> Complete Scan	<b>Digital Weighing ID:</b>
<b>Measurement Unit:</b> micron	

Size class for particle counting (ISO 16232-10)

Size X (µm)	Size Class	Metallic	Non-Metallic	Fibers	Total
5-15	B	123	19796	0	
15-25	C	106	2675	0	
25-50	D	11	1789	0	
50-100	E	15	1133	0	
100-150	F	2	247	56	
150-200	G	0	106	28	
200-400	H	0	77	24	
400-600	I	0	0	6	
600-1000	J	0	0	0	
1000-inf	K	0	0	0	
<b>Total Particle count</b>		<b>257</b>	<b>25823</b>	<b>114</b>	<b>26194</b>

Number of particles per 1000cm <sup>2</sup> or per 100 cm <sup>3</sup> more than	Up to and	Contamination level
0	0	00
0	1	0
1	2	2
2	3	3
4	8	3
8	10	4
16	32	5
32	64	6
64	130	7
130	250	8
250	500	9
500	1x10 <sup>3</sup>	10
1x10 <sup>3</sup>	2x10 <sup>3</sup>	11
2x10 <sup>3</sup>	4x10 <sup>3</sup>	12
4x10 <sup>3</sup>	8x10 <sup>3</sup>	13
8x10 <sup>3</sup>	16x10 <sup>3</sup>	14
16x10 <sup>3</sup>	32x10 <sup>3</sup>	15
32x10 <sup>3</sup>	64x10 <sup>3</sup>	16
64x10 <sup>3</sup>	130x10 <sup>3</sup>	17
130x10 <sup>3</sup>	250x10 <sup>3</sup>	18
250x10 <sup>3</sup>	500x10 <sup>3</sup>	19
500x10 <sup>3</sup>	1x10 <sup>6</sup>	20
1x10 <sup>6</sup>	2x10 <sup>6</sup>	21
2x10 <sup>6</sup>	4x10 <sup>6</sup>	22
4x10 <sup>6</sup>	8x10 <sup>6</sup>	23
8x10 <sup>6</sup>	16x10 <sup>6</sup>	24



<b>Contamination level</b>	<b>15</b>
----------------------------	-----------

<b>Max Metallic</b> 113.07 (L) X 109.49 (B)	<b>Max Non Metallic</b> 399.30 (L) X 299.81 (B)	<b>Max Fibre</b> 511.87 (L) X 50.91 (B)
--	--	--

<b>Analyzed area mm<sup>2</sup></b>	1385
<b>Total field area mm<sup>2</sup></b>	1711
<b>ISO 16232-7 (n) code</b>	B123/C106/D117/E15/F2/G0/H0/I0/J0/K0 B19796/C2675/D1789/E1133/F247/G106/H77/I0/J0/K0
<b>Result :</b>	

MILLIPORE TEST OBSERVATION			
Specification (mg)	Paper Weight (Before) (mg)	Paper Weight (After) (mg)	Milipore Value (mg)
<b>QA Remark:</b>			
<b>Inspected By:</b>		<b>Verified By:</b>	